

Level Plus[®]

Magnetostrictive Liquid-Level Transmitters with Temposonics[®] Technology

CIP Tanks Application Note

INDUSTRY: Pharmaceutical & biotechnology

PRODUCT: Cleaning agents



APPLICATION OVERVIEW

Automatic tank gauging of the level in the CIP tank where the concentrated cleaning agents are diluted to be used in the CIP process of the bioreactor and associated equipment.

KEY REQUIREMENTS

- High performance during agitation
- Sanitary (hygienic) wetted parts
- Accurate level measurement of product for dosing

MTS SENSORS SOLUTION SoClean Transmitter

- 3-A certified design to standard 74-06 sensors and sensor fittings and connections
- The SoClean model is not affected by agitation as the rigid 5/8" OD pipe is strong enough not to bend or vibrate during agitation
- Accurate level measurement of the product level to ± 1 mm (0.039 in.)



SoCLEAN[®]
LPS1C1C24F0XXXXU05000S

MTS Systems Corporation
Sensors Division
3001 Sheldon Drive
Cary, N.C. 27513
Phone: +1 919 677-0100
E-mail: info.us@mtssensors.com

Document Part Number: 550672 Revision E 11/2019

MTS, Temposonics and Level Plus are registered trademarks of MTS Systems Corporation in the United States; MTS SENSORS and the MTS SENSORS logo are trademarks of MTS Systems Corporation within the United States. These trademarks may be protected in other countries. All other trademarks are the property of their respective owners. Copyright © 2019 MTS Systems Corporation. No license of any intellectual property rights is granted. MTS reserves the right to change the information within this document, change product designs, or withdraw products from availability for purchase without notice. Typographic and graphics errors or omissions are unintentional and subject to correction. Visit www.mtssensors.com for the latest product information.